
**FRANÇAIS**

**NOTICE D'INSTRUCTIONS  
PROTHESE TOTALE DE CHEVILLE SALTO TALARIS**

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
**ENGLISH**

**INSTRUCTIONS FOR USE  
SALTO TALARIS TOTAL ANKLE PROSTHESIS**

**IMPORTANT:** The manufacturer recommends that all personnel responsible for handling and implanting the devices read and understand this information before use. The implantation of a joint prosthesis and its associated implants requires knowledge of anatomy, biomechanics and reconstructive surgery of the musculo-skeletal system and may be performed only by a qualified surgeon. The surgeon must operate in accordance with current information on the state of scientific progress and the art of surgery. Surgical training is required as a prerequisite to implant the Salto Talaris Total Ankle Prosthesis. The patient must be properly informed about the device and the information contained in the present instructions for use. **Caution:** The Federal (United States) Law restricts this device to sale, distribution and use by or on the order of a physician.

**1. Description:** The Salto Talaris Total Ankle is a semi-constrained prosthesis consisting of two mating components: a metal tibial base in association with a polyethylene insert and a metal talar resurfacing component. Accessories are provided to ensure proper implantation: Templates for each implant size, to be positioned on the X-rays to select the size of the implant prior to the surgery. Ancillary instruments are also provided:  
 • trial pieces for testing implantation during the surgery,  
 • instruments for the assembly and proper implanting of the prosthesis.  
 For a more detailed description of the implants and their utilization, please refer to the technical documentation, or contact your Tornier representative. It is essential to implant the Salto Talaris Total Ankle Prosthesis with the Tornier instrumentation specifically designed for this purpose. Tornier implants must be assembled using Tornier components defined as being compatible with one another. The selection of the appropriate implants can be made by using the recommendations of the surgical technique and the trial pieces and accessories supplied with the instrumentation. Symbols can be used to identify some implants (labeling or marking). They have the following meaning:  
**L** = left; **R** = right; **TH** = thickness.  
 Recommendations about how to associate components with each other are mentioned in the following drawing:

**2. Materials:** The constituent material of the Salto Talaris implants is labeled on the packaging. The tibial base and the talar component are manufactured from cobalt chromium alloy (CoCr) according to ISO standard 5832-4. The polyethylene insert is manufactured from implant grade ultra high molecular weight polyethylene (UHMWPE) according to ISO standard 5834-2. Tibial and talar implants are coated with plasma-spray porous titanium according to ASTM F1580.

**3. Intended use:** The Salto Talaris Total Ankle Prosthesis is intended for replacement of the ankle joint to reduce pain and restore ankle function compared with preoperative status.

**4. Indications for use:** The Salto Talaris Total Ankle Prosthesis is indicated as a total ankle replacement in primary or revision surgery for patients with ankle joints damaged by severe rheumatoid, post-traumatic, or degenerative arthritis. All components are intended for cemented use only.

- 5. Known contraindications to date:**
- Sepsis
  - Infection sequelae.
  - Systemic infection, fever and/or local inflammation.
  - Complete talar necrosis.
  - Insufficient quantity of bone stock or poor skin coverage around the ankle joint that would make the procedure unjustifiable.
  - Persisting skin lesion.
  - Important ligament laxity.
  - Severe osteoporosis.
  - Ankle arthrosis with malleolar exeresis.
  - Neuromuscular or mental disorders which might jeopardize fixation and postoperative care.
  - Neurobiologic diseases.
  - Nonfunctional lower limb muscles.
  - Complete loss of ankle collateral ligament.
  - Charcot's arthropathy.
  - Elevation of WBC count.
  - Distant foci of infection from genitourinary, pulmonary, skin and other sites, dental focus infection which may cause hematogenous spread to the implant site.
  - Bone immaturity.
  - Known allergy to one of the materials.
  - Patient pregnancy.

**6. Side-effects and possible complications:** The following are the most frequent adverse events after ankle arthroplasty:

- dislocation,
- infection,
- poor wound healing.

- loosening of components,
  - instability,
  - bone fracture,
  - secondary necrosis of the talus,
  - neuropathies,
  - disassembly or breakage of components,
  - possible metal sensitivity.
- 7. Warnings and cautions:**
- **Never re-use an implant, even if it seems to be in perfect condition, to prevent any risks of cross-contamination or a risk of reduced performances.**
  - **Never re-sterilize an implant delivered sterile.**
- The following conditions tend to adversely affect ankle replacement implants:
- Obesity or excessive patient weight;
  - Manual labor;
  - Active sports participation and/or high activity level;
  - Likelihood of falls;
  - Alcohol and/or drug addiction;
  - Other disabilities, as appropriate;
  - Poor bone stock;
  - Metabolic disorders or systemic pharmacological treatments leading to progressive deterioration of solid bone support for the implant (e.g. diabetes, steroid usage, immunosuppressive treatments);
  - Compromise of the ligaments or other supporting soft tissue structures such that they cannot withstand expected loads following arthroplasty, due to, for example, rheumatoid arthritis or other diseases affecting the quality of the soft tissue;
  - Severe deformities of the joint;
  - Tumors of the supporting bone structures;
  - Sensitivity, allergy or other reactions to implant materials;
  - Elevation of sedimentation rate unexplained by rheumatoid arthritis may adversely affect ankle replacement implants.
  - Inability of the patient to follow the surgeon's recommendations and the physical therapy program.
  - Pre-operatively: The surgeon must be fully conversant with all aspects of the surgical technique and know the indications and contraindications for this type of implant. The surgeon must have acquainted himself before the operation with the specific operative technique of the product which is available from the Tornier representative. As part of the preoperative examination, the surgeon must check that no biological, biomechanical or other factors may affect the correct conduct of the operation and the postoperative period. He must also check that the quality of the bone is satisfactory enough to support the implantation. An appropriate range of sizes must be available at the time of the

operation.

- Intra-operatively: The correct selection of the type and size of the implant appropriate to the patient and the positioning of the implant are extremely important. Never use an insert of a larger size than the size of the tibial plate. Check the proper anterior/posterior positioning of the tibial implant before impaction and that the insert lateralization is respected according to the side operated on. The use of trial pieces allows for the proper size selection of the implants. Frequent radioscopic checks allow the position of the prosthesis to be checked. The prostheses must not be used if their functional surfaces have been damaged or have undergone shock, abrasion, or other deterioration.
- Post-operatively: The surgeon must inform patients about:
  - precautions to take in daily life to guarantee maximum implant survival,
  - the fact that their weight and level of activity can affect the life span of the prosthesis.
 The Salto Talaris prosthesis has not been evaluated for safety and compatibility in the MR environment. The Salto Talaris prosthesis has not been tested for heating or migration in the MR environment. It is recommended that a regular postoperative follow-up is undertaken to detect early signs of wear, loosening of the prosthesis, etc., and to consider an appropriate course of action. Normal wear of the implant in respect of the state of knowledge at the time of its design cannot in any way be considered to constitute a malfunction or a deterioration in the characteristics of the implant. A suitable rehabilitation program must be designed and implemented specific to the patient.

**8. Storage and handling:** Implants must be stored in their original sealed packaging. The storage place must be away from humidity. Implants must not be exposed to direct sunlight, ionising radiation, extreme temperatures nor particular contamination. Implants must be handled with care to preserve integrity of the packaging.

**9. Packaging and sterilization:** The implants are supplied sterile (gamma radiation). An implant whose packaging is open or damaged or whose expiration date has passed must not be used. Every precaution must be taken to ensure sterility when opening the packaging of the implant and during implantation. Ancillary instruments may be supplied sterile. For handling and sterilization of non-sterile ancillary instruments, refer to the ancillary instruments instructions. The templates are supplied non-sterile and should not be sterilized.

For any other information regarding the ancillary instruments, refer to the instructions provided for this purpose.

**Interpretation of terms and symbols**

- ☒ Do not reuse
- ☒ Use by
- STERILE Irradiation
- LOT Lot number (sterilization)
- SN Serial number (manufacturing)
- REF Catalogue number
- ▲ See Instructions for use
- 📖 See technical documentation and/or surgical technique
- ⚠ Non sterile
- ☒ Do not use if package is damaged
- ☒ Do not re-sterilize
- ☒ Latex Free

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*The French text is the reference text. CE marking is only valid if it is also mentioned on the external package labeling*

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TORNIER, Inc